# Lucitone<sup>®</sup> *Fas-Por*<sup>™</sup> +

Pourable Denture Base

Plaque base coulable

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Prothesenkunststoff

Resina colable

**DIRECTIONS FOR USE**

**MODE D'EMPLOI**

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**Indications for use:** Lucitone® Fas-Por™+ Pourable Denture Base is an auto-polymerizing powder and liquid denture base material based on polymethyl methacrylate. This product is intended for the fabrication of complete or partial removable dentures or other removable appliances. The Lucitone Fas-Por+ Pourable Denture Base can also be used for repairs and relines.

### Contraindications:

1. Lucitone Fas-Por+ Pourable Denture Base is contraindicated for patients and users with a history of allergic reaction to methyl methacrylate monomer.

### Warnings:

1. Lucitone Fas-Por+ Pourable Denture Base contains polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wash thoroughly with soap and water after contact. If dermatitis or other symptoms persist, seek medical advice.
2. Avoid inhalation or ingestion. High vapor concentration can induce headache, nausea, and irritation of eyes and respiratory system. Liquid contact with eyes may cause possible corneal damage. Excessive long-term exposure may be associated with other more serious health effects. Monitor air quality per OSHA standards.

**Inhalation:** Move subject to fresh air. Give oxygen or artificial respiration as required.

**Ingestion:** Contact your regional Poison Control Center immediately.

**Eye Contact:** Flush eyes promptly with copious amounts of water for 15 minutes, and consult a physician. Wash skin with soap and water.

### Precautions:

1. Users with special skin problems, cuts, or abrasions should wear protective gloves.
2. When grinding prosthodontic resins, proper ventilation, masks, and vacuum systems should be used.
3. Store at 60°F-80°F (16°C-27°C), away from moisture and direct sunlight. The liquid contains unsaturated monomers which may polymerize prematurely if stored at excessively high temperature and/or in sunlight. The powder contains an organic polymerization initiator which may de-

grade if stored at excessively high temperature.

4. Lucitone Fas-Por+ Pourable Denture Base Liquid and Lucitone Fas-Por+ Pourable Denture Base Uncured resin are hazardous materials. Dispose of in accordance with Federal, State, and Local regulations.
5. Lucitone Fas-Por+ Pourable Denture Base Liquid contains methyl methacrylate monomer, a flammable liquid with a flash point of 50°F (10°C). Keep away from heat, sparks, and open flame.
6. Use liquid in a well ventilated area. Replace cap when not in use.
7. Do not use detergents to clear residual wax from casts. Detergent residue may impact the color uniformity.
8. Lucitone Fas-Por+ Pourable Denture Base resin's setup time is quick (approximately 1 minute). **After mixing powder and liquid, IMMEDIATELY pour into flask. Do not pour the denture resin after the material has thickened.**
9. **DO NOT ALLOW** Lucitone Fas-Por+ Pourable Denture Base to bench set prior to curing. Place into the pressure pot **IMMEDIATELY** after pouring.
10. Do not use \*Cetylcide-G® in solution to soak models.

### Adverse Reactions:

1. Corneal damage, headache, nausea, vomiting may occur with exposure to methyl methacrylate monomer. (See Warnings and Precautions statements.)
2. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals. Residual monomer in fully cured materials can be minimized by soaking the cured prosthesis in warm water for several days.
3. Particulates will be generated when grinding acrylate resins. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

---

**STEP-BY-STEP INSTRUCTIONS: PROPERTIES AND RECOMMENDED PROCEDURES** - The following table describes typical properties of Lucitone Fas-Por+ Pourable Denture Base, plus recommended temperatures and other important processing information. Further detail of this information is found later in the text.

Maximum Residual MMA Content .....	2.9% (% mass fraction)
Type and Class (according to ADA Spec. 12) .....	Type II, Class II
Storage temperature for powder and liquid .....	60°F-80°F (16°C-27°C)
Powder / Liquid ratio (Use Lucitone Fas-Por+ Liquid) .....	20 g (29 cc) / 14ml
Mixing time (time required to wet all particles) .....	15 seconds
Time to reach pourable stage @ 70°F ± 2°F .....	pour immediately
Pour Time .....	1 minute
Material used to prepare mold .....	hydrocolloid
Temperature of hydrocolloid when pouring .....	60°F-80°F (16°C-27°C)
Cure time, temperature, and pressure .....	pressure pot, 15 minutes @ 20-30psi, 106°F ± 2°F (41°C ± 1°C)

## FABRICATION OF FULL AND PARTIAL DENTURES

### HYDROCOLLOID INVESTMENT AND SPRUING

- Incorporate a slight taper into base of cast to facilitate removal from the hydrocolloid. Taper the cast so the base is the widest. **Soak the waxed stone cast in slurry water (saturated calcium sulfate dihydrate solution) for 10-20 minutes prior to pouring the hydrocolloid mold.** Do not use \*Cetylclde-G® in the solution. Follow the manufacturers instructions for preparation and use of the selected hydrocolloid product and flask. Bore sprue holes in the hydrocolloid investment as shown in Figure 1. **IMPORTANT:** When fabricating partials each acrylic area must have two sprues to avoid trapped air.

### WAX ELIMINATION AND SEPARATOR APPLICATION

- Boil off wax with clean boiling water. **DO NOT** use detergent to clean casts. Detergent residue may impact color uniformity. Also clean teeth with clean boiling water. Add mechanical retention to the teeth. Remove only enough tooth material to enhance retention. Excessive removal of tooth material will weaken the tooth and may lead to fracture.

Immediately after boil out, place the cast in a bowl of warm water for 5 minutes. Remove the cast from the water and allow to dry approximately 1-2 minutes.

Apply a thin coat of Al-Cote® Separating Agent to the cast and allow to dry until a slight sheen remains. Carefully reposition the teeth and cast back into the hydrocolloid mold. Reassemble the base and cover of the flask, being especially careful to avoid possible tooth and / or cast movement.

**MIXING** - To assure uniformity, shake powder jar (with a rocking motion) to prevent concentration of small particles at the bottom of jar. For the fabrication of 1 full denture **use 29 cc's (20 grams) of powder to 14 ml Lucitone Fas-Por+ Liquid.** **Do not use DENTSPLY Repair Material/ Lucitone Fas-Por E pourable denture base liquid.** Stir sufficiently (15 seconds) to assure wetting of all powder particles. Do not over-mix. **Pour immediately!**


**POURING** - Lucitone Fas-Por+ Pourable Denture Base should be **poured immediately after mixing.** Pour into one of the sprue openings with the flask positioned at a slight angle.

**IMPORTANT:** Pouring into more than one sprue opening may cause trapped air.

**POUR TIME** - The pour time of Lucitone Fas-Por+ Pourable Denture Base is approximately 1 minute. If the pour is done quickly this time enables the fabrication of 2 full dentures with one mix of material.

**IMPORTANT:** Immediately after pouring the dentures they must be placed into a pressure pot. **DO NOT BENCH SET.**

**CURING** - Check the water bath temperature in the pressure pot with an accurate thermometer. Place the flask in the water at 106°F ± 2°F (41°C ± 1°C) at 20-30 psi for 15 minutes. The water level should be at sprue height. **DO NOT FULLY SUBMERGE THE FLASK IN WATER.** If the water temperature is not maintained at approximately 106°F (41°C), the cure time will have to be increased to achieve full properties. For example, cure time should be 30 minutes at 73°F (23°C).

**Figure 1. Spruing Techniques**

**FINISHING AND POLISHING** - Lucitone Fas-Por+ Pourable Denture Base should be finished and polished using normal procedures.

## REPAIRING DENTURES

Lucitone Fas-Por+ Pourable Denture Base is recommended for repairing acrylic denture materials.

### PREPARATION OF REPAIR AREAS

Prepare the area by grinding a bevel on the edges. Wet the repair area with the Lucitone Fas-Por+ Liquid to enhance bonding.

### PROCESSING

When repairing acrylics Lucitone Fas-Por+ Pourable Denture Base should be cured in a pressure curing unit. The result is maximum density in the repaired area. Porosity and/or surface blanching may occur

if a pressure curing unit is not used. Submerge the repaired appliance (and cast as appropriate) in a warm water bath at 106°F (41°C ± 1°C) and 20-30 psi pressure for 15 minutes.

## INDIRECT RELINE PROCEDURES

Lucitone Fas-Por+ Pourable Denture Base is also recommended for relining poorly fitting dentures by the indirect method. A reline jig is recommended. Grind the tissue surface to relieve the denture. Soak the cast in slurry water for 10-20 minutes to eliminate air from the stone (prevents air bubbles in the acrylic). Wet the tissue side of the denture with the Lucitone Fas-Por+ Pourable Denture Base Liquid, then reseat the cast / denture on the jig. Immediately mix the Powder / Liquid (29cc / 14ml) and allow to bench set to reach the proper viscosity for the reline procedure. Process same as for repairs and dentures. A pressure pot is required to achieve the highest possible density.

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\*Cyticide-G® is not trademark of DENTSPLY
© 2009 DENTSPLY International Inc. All rights reserved.
**DENTSPLY International Inc.**
**DeguDent GmbH**